

# Installing Your Connect America Medical Alert System

The Connect America medical response system is an easy, self-install **medical alert system** that plugs into your phone just like an answering machine! Save time and money by self installing the system yourself, or have a friend, neighbor or relative do it for you. We simply ship you the unit and you follow the steps below to install and activate.

## Follow these 3 simple steps to self-install:

---



### STEP 1:

Place the Connect America **medical alert system** (the main communicator unit) in a central place of the home by the main phone and near a power outlet. In a home with multiple levels, we suggest the main floor and in an area away from high humidity. Plug the AC adapter into the power outlet.



### STEP 2:

Unplug your phone from the wall jack. Plug the end of your phone line into the back of the main Connect America communicator unit into the input marked "Phone". Using the phone cord supplied, plug the larger end into the input on the back of the main Connect America communicator unit marked "Line". Finally, plug the smaller end into the existing wall jack. Flip the switch on the back of the Connect America alert system to "On".



### STEP 3:

Call **1-800-325-1787** to activate your Connect America medical alert system. We will walk you through the activation procedure. When instructed to do so, you will hang up the phone and press your Personal Help button. Our Emergency Response Operator will greet you through the Connect America **medical alert system** to confirm your setup is complete!